

WONDFO BIOTECH
WeAreWorkingForYour **Health**

H.pylori Ag

FinecareTM

H.pylori Antigen Test

**Rapid and accurate indication of
H.pylori infection**

Wondfo

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City,
Luogang District, 510663 Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com.cn
E-mail: sales@wondfo.com.cn



Rev. A2 Rel.: 2023/01/29



H.pylori Antigen Test

- ✓ Aid in the diagnosis of H.pylori infection
- ✓ Qualitative determination of H.pylori Ag in human feces
- ✓ Good clinical sensitivity and specificity
- ✓ Room temperature storage up to 24 months
- ✓ User-friendly, non-invasive and fast detection method
- ✓ Stool sample supported, perfect for the elderly, non-compliant children and pregnant women.
- ✓ Convenient all-in-one sampler included

CLINICAL APPLICATION

Diagnosis

Aid in the diagnosis of H.pylori infection in adult patients with symptoms of gastrointestinal disorders.

Monitoring

ACG: Whenever H.pylori infection is identified and treated, testing to prove eradication should be performed using a urea breath test, fecal antigen test or biopsy-based testing at least 4 weeks after the completion of antibiotic therapy.

«ACG Clinical Guideline: Treatment of Helicobacter pylori Infection»
*ACG: American College of Gastroenterology

Application scenarios



GP Office



Central Laboratory



Clinics



Pharmacy



Medical Examination Center




Gastroenterology



ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Reaction Time	Storage	Sensitivity	Specificity	CV
H.pylori Ag Test	W270P0001	Human feces	10 mins	4-30 °C	99.0%	98.8%	≤10%

TEST PROCEDURE



ADAPTATION INSTRUMENTS

1. Finecare™ FIA Meter Plus (FS-113)
2. Finecare™ II Plus SE (FS-114)
3. Finecare™ III Plus (FS-205)

