

WONDFO BIOTECH
WeAreWorkingForYour **Health**

Wondfo

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City,
Luogang District, 510663 Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com.cn
E-mail: sales@wondfo.com.cn



Rev. A1 Rel.: 2022/03/28

ASO

Anti-streptolysin O


Finecare™

ASO Rapid Quantitative Test

For diagnose rheumatic fever and glomerulonephritis
induced by streptococci infection

ASO Rapid Quantitative Test

- ✓ Whole blood samples supported
- ✓ Results ready in 5 minutes
- ✓ Room temperature storage up to 24 months



CLINICAL APPLICATION

To help determine whether you have had a recent strep infection with the bacteria group A Streptococcus;


To help reflect the infection status and extent of the disease;

To help diagnose complications resulting from a strep infection such as rheumatic fever or glomerulonephritis

ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Loading Volume	Reaction Time	Detection Range	Storage
ASO	W297P0001	whole blood, serum and plasma	10 µL	75 µL	5 minutes	20~1000 IU/mL	4~30°C for 24 months

TEST PROCEDURE



ADAPTATION INSTRUMENTS

1. Finecare™ FIA Meter Plus (FS-113)
2. Finecare™ II Plus SE (FS-114)
3. Finecare™ III Plus (FS-205)

