



WONDFO BIOTECH
WeAreWorkingForYour **Health**

IgE

Immunoglobulin E

Finecare™

Total IgE Rapid Quantitative Test

Helping to diagnose allergic diseases



Guangzhou Wondfo Biotech Co., Ltd.
No.8 Lizhishan Road, Science City, Luogang District,
510663, Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com
E-mail: sales@wondfo.com.cn



Rev. A1 Rel.: 2022/03/28

Total IgE Rapid Quantitative Test

- Whole blood samples supported
- Instant results in 10 minutes
- Room temperature storage up to 24 months




CLINICAL APPLICATION

- ✓ Helping to diagnose allergic diseases
- ✓ Predicting the development of allergic disorders
- ✓ Auxiliary diagnosis of allergic rhinitis
- ✓ Auxiliary diagnosis of asthma

ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Loading Volume	Reaction Time	Detection Range	Storage
Total IgE	W298 P0001	whole blood, serum and plasma	75 μ L	75 μ L	10 minutes	1~1000 IU/mL	4~30°C for 24 months

TEST PROCEDURE



ADAPTATION INSTRUMENTS

1. Finecare™ FIA Meter Plus (FS-113)
2. Finecare™ II Plus SE (FS-114)
3. Finecare™ III Plus (FS-205)

