

WONDFO BIOTECH
WeAreWorkingForYour **Health**

Wondfo

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City,
Luogang District, 510663 Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com.cn
E-mail: sales@wondfo.com.cn



Rev. A1 Rel.: 2022/01/11

Wondfo reserves all rights. All contents of this brochure are copyrighted by Wondfo.

Wondfo



25-OH-VD
more than just a vitamin **D**

FinecareTM

VITAMIN D RAPID QUANTITATIVE TEST

Only four steps, easy to use !

Finecare™


Vitamin D Rapid Quantitative Test

- ▶ An effected biomarker for screening and diagnosing vitamin D deficiency
- ▶ Room temperature storage up to 24 months
- ▶ Results in less than 25 minutes
- ▶ Simple four steps, Easy to use
- ▶ Minimum training required



CLINICAL APPLICATION

Diagnosing Vitamin D deficiency
 Diagnosing childhood rickets
 Diagnosing adult osteomalacia
 Diagnosing hypercalcemia




CLINICAL REFERENCE

< 20 ng/mL Deiciency ;
 20~30 ng/mL Insuiciency
 30~100 ng/mL Suiciency ;
 > 100 ng/mL Toxicity

ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Loading Volume	Reaction Time	Detertion Range	Storage
Vitamin D	W241P0002	Serum / Plasma	75 µL	75 µL	25 minutes	5~100 ng/mL	4~30 C for 24 months

TEST PROCEDURE



ADAPTATION INSTRUMENTS

1. Finecare™ FIA Meter Plus (FS-113)
2. Finecare™ II Plus SE (FS-114)
3. Finecare™ III Plus (FS-205)



Point of Care (GP, Doctor's Office, Clinics)



Laboratory



Pediatric Department



Orthopedics



Maternal and Child Health



Geriatrics