



H. Pylori

One Step H. Pylori Serum/Plasma Test

Wondfo One Step H. Pylori Serum/Plasma Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies to H. pylori in serum or plasma to aid in the diagnosis of H. pylori infection.

Specimen: Serum/Plasma

Detection Target: H. Pylori antibody

Detection Time: 15mins

Sensitivity: 98.35%

Specificity: 98.23%



Guangzhou Wondfo Biotech Co., Ltd.

No.8 Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R.China

Tel: (+86)400-830-8768 E-mail: sales@wondfo.com.cn

Website: en.wondfo.com.cn


Rev. A1 Rel.:2022/06/30

WONDFO BIOTECH


WeAreWorkingForYourHealth

Test Procedure


Step 1



Step 2



Step 3



Performance Characteristics

Wondfo One Step H. Pylori Serum/Plasma Test	Commercial kit		Total
	Positive	Negative	
Positive	476	10	486
Negative	8	556	564
Total	484	566	1050

Sensitivity of Wondfo One Step H.Pylori Serum/Plasma Test: 98.35%

Specificity of Wondfo One Step H.Pylori Serum/Plasma Test: 98.23%

Order Information

Catalog No.	Product Name	Package Size	Certification
W38-C	One Step H.Pylori Serum/Plasma Test	25T/kit	CE
W038P0043	One Step H.Pylori Serum/Plasma Test	40T/kit	CE