

Wondfo

HIV-HCV-HBsAg-Syphilis Combo Test



Wondfo One Step HIV-HCV-HBsAg-Syphilis Combo Whole Blood/Plasma/Serum Test is a rapid immunochromatographic test for visual detection of HIV1/2 antibodies, Hepatitis C Virus (HCV) antibody, hepatitis B surface antigen (HBsAg) and Treponema Pallidum (TP) antibody in human whole blood/plasma/serum.



Specimen | Whole blood/serum/plasma

Time to results | 15 mins

Storage temperature | 4~30°C

Shelf-life | 24 months

Guangzhou Wondfo Biotech Co., Ltd.

No.8 Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R.China




Tel: (+86)400-830-8768 E-mail: sales@wondfo.com.cn Website: en.wondfo.com.cn

Rev. A1 Rel.:2022/06/08





01. TEST PROCEDURE IN 3 SIMPLE STEPS

For whole blood sample testing



For serum/plasma sample testing



02. PERFORMANCE CHARACTERISTICS

Items	Sensitivity	Specificity
HIV	99.70%	100%
HCV	98.71%	99.33%
HBsAg	98.97%	99.78%
Syphilis	98.47%	99.68%

03. ORDER INFORMATION

Product name	Catalog no.	Pack size
One Step HIV-HCV-HBsAg-Syphilis Combo Whole Blood/Plasma/Serum Test	W026-C	25T/kit