

### DOA Urine Tests

Qualification: **CE / FDA / NMPA**


Various formats: **Cassette, Panel, T-Cup, Q-Cup**

Features: **Optional adulteration test**

Storage TEM: **4°C~30°C**

Shelf life: **24 months**

Principle: **Competitive immunoassay**



#### • One Step DOA Single Tests (22 items are available)

Item	Sample	Cutoff (ng/mL)
Amphetamine (AMP)	Urine	1000
Barbiturates (BAR)	Urine	300
Benzodiazepines (BZO)	Urine	300
Cocaine (COC)	Urine	300 / 150
Methamphetamine (mAMP/MET)	Urine	1000
Methylenedioxymethamphetamine; Ecstasy (MDMA)	Urine	500
Buprenorphine (BUP)	Urine	10
Morphine (MOP)	Urine	300
Methadone (MTD)	Urine	300
Opiate (OPI)	Urine	2000
Phencyclidine (PCP)	Urine	25

Item	Sample	Cutoff (ng/mL)
Tricyclic Antidepressants (TCA)	Urine	1000
Marijuana (THC)	Urine	50 / 25
Oxycodone (OXY)	Urine	100
Ketamine (KET)	Urine	1000
Propoxyphene (PPX)	Urine	300
Ethylenediamine-dimethyphosphinic (EDDP)	Urine	100
Tramadol (TRA)	Urine	200
Fentanyl (FTY)	Urine	20
Synthetic Cannabis (K2)	Urine	50
Cotinine (COT)	Urine	100
Ethyl Glucuronide (ETG)	Urine	500

#### • One Step DOA Multi-Tests (Various combinations are available)

Product	Test items	Format	Specimen	Sensitivity	Qualification
DOA Multi-Test Panel	2-12 Items; 20+ for options	Panel	Urine	Variable	CE FDA
DOA Multi-Test Panel	13-16 Items; 20+ for options	Panel	Urine	Variable	CE FDA
DOA Multi-Test Cup	2-18 Items; 20+ for options	T-Cup	Urine	Variable	CE FDA
DOA Multi-Test Cup	2-18 Items; 20+ for options	Q-Cup	Urine	Variable	CE



Cup



Panel

## DOA Oral Fluid Tests

Principle: **Colloidal gold immunochromatography**

Feature: **Compact size with up to 2-10 items in T-Cube**

Advantage: **High sensitivity**

Storage TEM: **4°C~30°C**

Shelf life: **24 months**



Cat.No	Product	Specimen	Format	Qualification
W602~610-CU6	<b>Multi-Drug Oral Fluid</b>	Oral Fluid	T-Cube	CE

### • Multi-DOA T-Cube Test Items

Item	Sample	Cutoff (ng/mL)	Item	Sample	Cutoff (ng/mL)
<b>Amphetamine (AMP)</b>	Oral Fluid	50	<b>Methamphetamine (MET)</b>	Oral Fluid	50
<b>Barbiturate (BAR)</b>	Oral Fluid	60	<b>Methadone (MTD)</b>	Oral Fluid	30
<b>Buprenorphine (BUP)</b>	Oral Fluid	5	<b>Morphine (MOP)</b>	Oral Fluid	40
<b>Benzodiazepines (BZO)</b>	Oral Fluid	30	<b>Opiates (OPI)</b>	Oral Fluid	40
<b>Cocaine (COC)</b>	Oral Fluid	20	<b>Phencyclidine (PCP)</b>	Oral Fluid	10
<b>Synthetic Marijuana; Synthetic Cannabis (K2)</b>	Oral Fluid	5	<b>Marijuana; Tetrahydrocannabinol; Cannabis (THC)</b>	Oral Fluid	40
<b>Methylenedioxymethamphetamine; Ecstasy (MDMA)</b>	Oral Fluid	100	<b>Oxycodone(OXY)</b>	Oral Fluid	20

### • Other Rapid Test (Various combinations are available)

Cat.No	Product	Specimen	Format	Qualification
W53-S	<b>Alcohol Saliva Test</b>	Oral Fluid	Strip	CE FDA

Guangzhou Wondfo Biotech Co., Ltd.

No. 8 Lizhishan Road, Science City, Huangpu District, 510700, Guangzhou, P.R. China

E-mail: sales@wondfo.com.cn

Tel: (+86) 400-830-8768

Website: en.wondfo.com

Wondfo reserves all rights. All contents of this brochure are copyrighted by Wondfo.

