



WONDFO BIOTECH  
WeAreWorkingForYour **Health**

# Fertility

## Finecare™

### Fertility Rapid Quantitative Test

one-stop solution for sex hormones testing !



**Guangzhou Wondfo Biotech Co., Ltd.**  
No. 8 Lizhishan Road, Science City,  
Luogang District, 510663 Guangzhou, P.R.China  
Tel: (+86) 400-830-8768  
Website: en.wondfo.com.cn  
E-mail: sales@wondfo.com.cn



Rev. A1 Rel.: 2022/08/25

# Finicare™

## Fertility Rapid Quantitative Test

- ✓ Direct detection of sex hormone levels
- ✓ Help diagnosis PCOS, AUB, Infertility
- ✓ Eight sex hormone items, one-stop solution
- ✓ Room temperature storage up to 24 months

### CLINICAL APPLICATION

- **Children:** auxiliary diagnosis of precocious puberty;
- **Adolescents:** auxiliary diagnosis of sexual developmental delay;
- **Reproductive period:** auxiliary diagnosis of infertility, prenatal and postnatal care;
- **Menopause:** auxiliary diagnosis of changes in hormone levels in women during menopause

### APPLICATION SCENARIOS



### ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Loading Volume	Detection Range	Reaction Time	Storage
LH	W244	Whole blood/ Serum/Plasma	75 µL	75 µL	1~100 mIU/mL	15 minutes	4~30°C for 24 months
FSH	W245		75 µL	75 µL	1~100 mIU/mL		
Prog	W233		75 µL	75 µL	1.5~60 ng/mL		
PRL	W246P0001		75 µL	75 µL	1~200 ng/mL		
E2	W247-C4P1-E		75 µL	75 µL	9~3000 pg/mL		
Testosterone	W248-C4P1-E		75 µL	75 µL	0.2~15 ng/mL		
AHM	W257P0001		75 µL	75 µL	0.1~16 ng/mL		
β-hCG	W225		75 µL	20 µL	2~200,000 mIU/mL		

### TEST PROCEDURE

**Report**

Age : 27 years  
Gender: Female

**Result Display**

Results :  
AMH 0.2 ng/mL ↓

**Normal Reference Interval:**

Gender	Age (years old)	Reference interval (ng/mL)
Adult Male	-	0.92-13.89
Female	20-29	0.88-10.35
	30-39	0.31-7.86
	40-50	≤ 5.07

### ADAPTATION INSTRUMENTS

1. Finicare™ FIA Meter Plus (FS-113)
2. Finicare™ II Plus SE (FS-114)
3. Finicare™ III Plus (FS-205)

